

# Analysenzertifikat

## Certificate of Analysis

### SERION KBR Antigen / CFT Antigen VARICELLA-ZOSTER VIRUS (VZV)

Best.-Nr. / Order No.	1191
Chargen-Nr. / Lot No.	KK0052
IFU-Version / IFU version	21-19/07-1
Verfallsdatum / Expiry Date	2029-11-30
Antigen / Antigen	Varizella-Zoster-Virus (Ellen)
Kultursystem / Culture System	HEL299 Zellen / cells
Präparation / Preparation	infizierte Zellen und Zellkulturmedium (konzentriert und gereinigt) / infected cells and cell culture medium (concentrated and purified)
Inaktivierung / Inactivation	Physikalische Inaktivierung / physical inactivation
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:12



Hersteller / Manufacturer:  
 Institut Virion\Serion GmbH  
 Friedrich-Bergius-Ring 19  
 D-97076 Würzburg  
 Germany  
 Tel. +49 931 3045 0  
 Fax +49 931 3045 100  
 E-Mail: info@virion-serion.de  
 Internet: www.virion-serion.de