

Analysenzertifikat

Certificate of Analysis

SERION KBR Antigen / CFT Antigen PICORNA VIRUS POOL (EV, CAV, CBV)

Best.-Nr. / Order No.	1126
Chargen-Nr. / Lot No.	KL0037
IFU-Version / IFU version	21-19/07-1
Verfallsdatum / Expiry Date	2030-03-31
Antigen / Antigen	Coxsackie-Viren / Coxsackie viruses A9, B1, B2, B3, B4, B5, B6 (PB, Conn, Ohio 1, Nancy, JVB, Faulkner, Schmitt) Echo-Viren / Echo viruses 4, 6, 9, 14, 24, 30 (Pesascek, D' Amori, Hill, Tow, DeCamp, Bastianni)
Kultursystem / Culture System	LLC-MK2 Zellen / cells Vero-Zellen / cells
Präparation / Preparation	Infizierte Zellen und Zellkulturmedium (konzentriert und gereinigt) Infected cells and cell culture medium (concentrated and purified)
Inaktivierung / Inactivation	β -Propiolakton / UV-C β -Propiolactone / UV-C
Konservierung / Preservation	0,01 % Thiomersal / 0,01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:6



Hersteller / Manufacturer:
 Institut Virion\Serion GmbH
 Friedrich-Bergius-Ring 19
 D-97076 Würzburg
 Germany
 Tel. +49 931 3045 0
 Fax +49 931 3045 100
 E-Mail: info@virion-serion.de
 Internet: www.virion-serion.de