

Analysenzertifikat

Certificate of Analysis

SERION KBR Zusatzreagenz / CFT Supplementary Reagent
HÄMOLYTISCHER AMBOZEPTOR / HEMOLYTIC AMBOCEPTOR

Best.-Nr. / Order no.	9002
Chargen-Nr. / Lot no.	KL0042
IFU-Version / IFU version	21-19/07-1
Verfallsdatum / Expiry Date	2023-04-30
Material / Material	Kaninchen-anti-Hammel-Erythrozyten-Serum / rabbit anti-sheep erythrocyte serum
Lagerung / Storage	2 °C - 8 °C
Gebrauchsverdünnung / Working Dilution	1 : 2500



Hersteller / Manufacturer:
 Institut Virion\Serion GmbH
 Friedrich-Bergius-Ring 19
 D-97076 Würzburg
 Germany
 Tel. +49 931 3045 0
 Fax +49 931 3045 100
 E-Mail: info@virion-serion.de
 Internet: www.virion-serion.de