

# Analysenzertifikat

## Certificate of Analysis

### SERION KBR Antigen / CFT Antigen ECHOVIRUS POOL (4,6,9,14,24,30)

Best.-Nr. / Order No.	1180
Chargen-Nr. / Lot No.	KL0044
IFU-Version / IFU version	21-19/07-1
Verfallsdatum / Expiry Date	2030-05-31
Antigen / Antigen	Echoviren Typ / Echoviruses type 4, 6, 9, 14, 24, 30 (Stämme / strains: Pesascek, D'Amori, Hill, Tow, DeCamp, Bastianni)
Kultursystem / Culture System	Verozellen / Vero cells
Präparation / Preparation	infizierte Zellen und Zellkulturmedium (konzentriert und gereinigt) / infected cells and cell culture supernatant (concentrated and purified)
Inaktivierung / Inactivation	beta-Propiolacton und UV Behandlung / beta-propiolactone and UV treatment
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:6



Hersteller / Manufacturer:  
 Institut Virion\Serion GmbH  
 Friedrich-Bergius-Ring 19  
 D-97076 Würzburg  
 Germany  
 Tel. +49 931 3045 0  
 Fax +49 931 3045 100  
 E-Mail: info@virion-serion.de  
 Internet: www.virion-serion.de  
 FB-IV-047 D