

# Analysenzertifikat

## Certificate of Analysis

### SERION KBR Positivkontrolle / CFT Positive Control COXSACKIEVIRUS POOL (A9, B1-B6)

|   |                          |
|---|--------------------------|
| Best.-Nr. / Order No.                     | 3178                     |
| Chargen-Nr. / Lot No.                     | KL0066                   |
| IFU-Version / IFU version                 | 21-19/07-1               |
| Verfallsdatum / Expiry Date               | 2030-09-30               |
| Material / Material                       | Humanserum / human serum |
| Lagerung / Storage                        | 2 °C - 8 °C              |
| Rekonstitution / Reconstitution           | 0,1 ml Aqua dest.        |
| Gebrauchsverdünnung<br>/ Working Dilution | 1:10                     |
| KBR Titer / CFT Titre                     | 1:40                     |



Hersteller / Manufacturer:  
 Institut Virion\Serion GmbH  
 Friedrich-Bergius-Ring 19  
 D-97076 Würzburg  
 Germany  
 Tel. +49 931 3045 0  
 Fax +49 931 3045 100  
 E-Mail: info@virion-serion.de  
 Internet: www.virion-serion.de