

Specifications SERION easyCSF

PC

In order for SERION easyCSF to run on your PC, various hardware and software requirements must be met, which are listed below.

Operating system:	Win 10 Professional 64bit
Operating system language:	German or English (US)
Country-specific format:	German or English (US)
CPU:	Min. 2 x 1.8 GHz
RAM:	Min. 4 GB RAM

Operating system

SERION easyCSF can be used under Windows® 10. The application was validated under Windows 10, taking into account German and English (US) language and country (format) settings.

- Windows 10
- Microsoft® .NET Core 3.1.9 runtime or newer must be installed