

Qualitätskontrollzertifikat

Quality Control Certificate

SERION ELISA *classic* COXIELLA BURNETII PHASE 2 IgM

Best.-Nr. / Order No.	ESR1312M
Chargen-Nr. / Lot No.	EM0162
Verfallsdatum / Expiry Date	2023-06-30
IFU-Version / IFU-Version	NEW! 131-15

Kontrollseren / Control sera	Chargen-Nr. / Lot No	Gemessener Wert / Measured value OD 405 nm	Gültigkeitsbereich / Validity range
Positivkontrolle / Positive control	ECM0295	0,85	> 0,43
Cut off-Serum / Cut off serum	ECM0296	0,39	0,20 – 0,78

Negativkontrolle / Negative Control Chargen-Nr. / Lot No.	COXIELLA BURNETII PHASE 2 IgM ECM0294
Teststreifen / Strips Chargen-Nr. / Lot No.	COXIELLA BURNETII PHASE 2 IgM ECM0247
Konjugat / Conjugate Chargen-Nr. / Lot No.	IgM + KJM018

22.07.2021

Prüfdatum / Date of control

Hersteller / Manufacturer:
 Institut Virion\Serion GmbH
 Friedrich-Bergius-Ring 19
 D-97076 Würzburg
 Germany
 Tel. +49 931 3045 0
 Fax +49 931 3045 100
 E-Mail: info@virion-serion.de
 Internet: www.virion-serion.de