


Qualitätskontrollzertifikat

Quality Control Certificate

SERION ELISA *agile*
SARS-CoV-2 IgM



Best.-Nr. / Order No.	ESR400M
Chargen-Nr. / Lot No.	AM030
Verfallsdatum / Expiry Date	2022-12-31
IFU-Version / IFU-Version	NEW! a400M-2

Kontrollseren / Control sera	Chargen-Nr. / Lot No	Gemessener Wert / Measured value	Gültigkeitsbereich / Validity range
Positivkontrolle / Positive control	AM015	0,577	> 0,289
Cut off-Serum / Cut off serum	AM016	0,242	0,121-0,484

Negativkontrolle / Negative Control Chargen-Nr. / Lot No.	AM014
Teststreifen / Strips Chargen-Nr. / Lot No.	AM017
Konjugat / Conjugate Chargen-Nr. / Lot No.	AM019
RF Absorbent (Z23) / RF Absorbent (Z23) Chargen-Nr. / Lot No.	AL037

28.10.2021

Prüfdatum / Date of control

Hersteller / Manufacturer:
 Institut Virion\Serion GmbH
 Friedrich-Bergius-Ring 19
 D-97076 Würzburg
 Germany
 Tel. +49 931 3045 0
 Fax +49 931 3045 100
 E-Mail: info@virion-serion.de
 Internet: www.virion-serion.de

FB-IV-312B