

## Analysenzertifikat Certificate of Analysis

**SERION KBR Kontrollantigen / CFT Control Antigen  
MUMPS VIRUS/PAROTITIS/OREILLONS**

|  |   |
|--|---|
| Best.-Nr. / Order No.                  | 2125  |
| Chargen-Nr. / Lot No.                  | KN0003  |
| IFU-Version / IFU version              | 22-20/02-1  |
| Verfallsdatum / Expiry Date            | 2028-03-31  |
| Kultursystem / Culture System          | Verozellen / vero cells   |
| Präparation / Preparation              | Konzentrierter Zellkulturüberstand und alkalisch extrahierte Zellen /<br>Concentrated cell culture supernatant and cells extracted with alkaline buffer |
| Konservierung / Preservation           | 0,01 % Thiomersal / 0.01 % thiomersal   |
| Lagerung / Storage                     | 2 °C - 8 °C   |
| Rekonstitution / Reconstitution        | 1 ml Aqua dest.   |
| Gebrauchsverdünnung / Working Dilution | 1:10  |



Hersteller / Manufacturer:  
 Institut Virion\Serion GmbH  
 Friedrich-Bergius-Ring 19  
 D-97076 Würzburg  
 Germany  
 Tel. +49 931 3045 0  
 Fax +49 931 3045 100  
 E-Mail: info@virion-serion.de  
 Internet: www.virion-serion.de