

Qualitätskontrollzertifikat

Quality Control Certificate

SERION ELISA *classic* Adenovirus IgM

Best.-Nr. / Order No.	ESR128M
Chargen-Nr. / Lot No.	EL0185
Verfallsdatum / Expiry Date	2022-10-31
IFU-Version / IFU-Version	NEW! 128-16

Kontrollseren / Control sera	Chargen-Nr. / Lot No	Gemessener Wert / Measured value OD 405 nm	Gültigkeitsbereich / Validity range
Positivkontrolle / Positive control	ECL0443	1,13	> 0,57
Cut off-Serum / Cut off serum	ECL0454	0,23	0,12 – 0,46

Negativkontrolle / Negative Control	Adenovirus IgM
Chargen-Nr. / Lot No.	ECL0442
Teststreifen / Strips	Adenovirus IgM
Chargen-Nr. / Lot No.	ECL0437
Konjugat / Conjugate	IgM ++
Chargen-Nr. / Lot No.	KJL010

03.12.2020

Prüfdatum / Date of control

Hersteller / Manufacturer:
 Institut Virion\Serion GmbH
 Friedrich-Bergius-Ring 19
 D-97076 Würzburg
 Germany
 Tel. +49 931 3045 0
 Fax +49 931 3045 100
 E-Mail: info@virion-serion.de
 Internet: www.virion-serion.de