

## Analysenzertifikat Certificate of Analysis

### SERION KBR Antigen / CFT Antigen LEPTOSPIRA POMONA

Best.-Nr. / Order No.	9100
Chargen-Nr. / Lot No.	SKI.DS
Verfallsdatum / Expiry Date	2028-10
Antigen / Antigen	Leptospira interrogans, Serovar pomona
Kultursystem / Culture System	Flüssigkultur / broth culture
Präparation / Preparation	mit Ultraschall behandeltes Zellmaterial / sonicated cellular material
Inaktivierung / Inactivation	0,1 % Thiomersal / 0.1 % thiomersal
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:8



Hersteller / Manufacturer:  
 Institut Virion\Serion GmbH  
 Friedrich-Bergius-Ring 19  
 D-97076 Würzburg  
 Germany  
 Tel. +49 931 3045 0  
 Fax +49 931 3045 100  
 E-Mail: info@virion-serion.de  
 Internet: www.virion-serion.de