

# Qualitätskontrollzertifikat

## Quality Control Certificate

### SERION ELISA *classic* COXIELLA BURNETII PHASE 2 IgM

Best.-Nr. / Order No.	ESR1312M
Chargen-Nr. / Lot No.	EN0081
Verfallsdatum / Expiry Date	2024-03-31
IFU-Version / IFU-Version	131-15

Kontrollseren / Control sera	Chargen-Nr. / Lot No	Gemessener Wert / Measured value OD 405 nm	Gültigkeitsbereich / Validity range
Positivkontrolle / Positive control	ECN0154	0,80	> 0,40
Cut off-Serum / Cut off serum	ECN0155	0,39	0,20 – 0,78

Negativkontrolle / Negative Control Chargen-Nr. / Lot No.	COXIELLA BURNETII PHASE 2 IgM ECN0153
Teststreifen / Strips Chargen-Nr. / Lot No.	COXIELLA BURNETII PHASE 2 IgM ECM0247
Konjugat / Conjugate Chargen-Nr. / Lot No.	IgM + KJN035

22.04.2022

Prüfdatum / Date of control

Hersteller / Manufacturer:  
 Institut Virion\Serion GmbH  
 Friedrich-Bergius-Ring 19  
 D-97076 Würzburg  
 Germany  
 Tel. +49 931 3045 0  
 Fax +49 931 3045 100  
 E-Mail: info@virion-serion.de  
 Internet: www.virion-serion.de