

# Analysenzertifikat

## Certificate of Analysis

### SERION KBR Antigen / CFT Antigen TOXOPLASMA GONDII

Best.-Nr. / Order No.	1331
Chargen-Nr. / Lot No.	KN0014
IFU-Version / IFU version	22-20/02-1
Verfallsdatum / Expiry Date	2032-08-31
Antigen / Antigen	Toxoplasma gondii Stamm / strain RH Toxoplasma gondii Stamm / strain BK
Kultursystem / Culture System	HEL-299
Präparation / Preparation	geerntete und durch Zentrifugation partiell gereinigte Tachyzoiten in Puffer  harvested and partially purified tachyzoites in buffer
Inaktivierung / Inactivation	Hitzeinaktivierung / heat inactivation 30 min, 56°C
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:8



Hersteller / Manufacturer:  
 Institut Virion\Serion GmbH  
 Friedrich-Bergius-Ring 19  
 D-97076 Würzburg  
 Germany  
 Tel. +49 931 3045 0  
 Fax +49 931 3045 100  
 E-Mail: info@virion-serion.de  
 Internet: www.virion-serion.de  
 FB-IV-047 E