

Analysenzertifikat Certificate of Analysis

SERION KBR Kontrollantigen / CFT Control Antigen
MASERN VIRUS / MEASLES / ROUGEOLE

Best.-Nr. / Order No.	2190
Chargen-Nr. / Lot No.	KO0003
IFU-Version / IFU version	22-20/02-1
Verfallsdatum / Expiry Date	2028-03-31
Kultursystem / Culture System	Verozellen / vero cells
Präparation / Preparation	Konzentrierter Zellkulturüberstand und alkalisch extrahierte Zellen concentrated cell culture supernatant and cells extracted with alkaline buffer
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:10



Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: info@virion-serion.de
Internet: www.virion-serion.de