

Analysenzertifikat Certificate of Analysis

**SERION KBR Kontrollantigen / CFT Control Antigen
COXSACKIEVIRUS POOL (A9, B1-B6)**

Best.-Nr. / Order No.	2178
Chargen-Nr. / Lot No.	KO0007
IFU-Version / IFU version	22-20/02-1
Verfallsdatum / Expiry Date	2027-01-31
Kultursystem / Culture System	Vero Zellen / Vero cells und / and LLC-MK2 Zellen / LLC-MK2 cells
Präparation / Preparation	Konzentrierter Zellkulturüberstand und alkalisch extrahierte Zellen / Concentrated cell culture supernatant and cells extracted with alkaline buffer
Konservierung / Preservation	0,01 % Thiomersal / 0,01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:10



Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: info@virion-serion.de
Internet: www.virion-serion.de