

SAFETY INFORMATION

Voluntary safety information following the Safety Data Sheet format according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 26 Mar 2024

Print date: 26 Mar 2024

Version: 1.1

virion\serion

Page 1/8

SERION ELISA classic Kit

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Trade name/designation:

SERION ELISA classic Kit

Article No.:

ESR*

1.2. Relevant identified uses of the substance or mixture and uses advised against

Use of the substance/mixture:

The product is intended for professional use.

1.3. Details of the supplier of the safety data sheet

Supplier (manufacturer/importer/only representative/downstream user/distributor):

Institut Virion\Serion GmbH

Productmanagement

Friedrich-Bergius-Ring 19

97076 Würzburg

Germany

Telephone: +49 931 3045 0

Telefax: +49 931 3045 100

E-mail: product.safety@virion-serion.de

Website: www.virion-serion.de

E-mail (competent person): product.safety@virion-serion.de

* 1.4. Emergency telephone number

Institut Virion\Serion GmbH, Friedrich-Bergius-Ring 19; 97076 Würzburg, Germany, +49 931 3045 0 (Only available during office hours.)

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No 1272/2008 [CLP]

The mixture is classified as not hazardous according to regulation (EC) No 1272/2008 [CLP].

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

According to EC directives or the corresponding national regulations the product does not have to be labelled.

Hazard components for labelling:

2-methyl-2H-isothiazol-3-one; hydrogen chloride

Hazard statements: none

Supplemental hazard information

EUH208

Contains 2-methyl-2H-isothiazol-3-one. May produce an allergic reaction.

Precautionary statements: none

2.3. Other hazards

No data available

SECTION 3: Composition/information on ingredients

3.2. Mixtures

No data available

SAFETY INFORMATION

Voluntary safety information following the Safety Data Sheet format according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 26 Mar 2024

Print date: 26 Mar 2024

Version: 1.1

virion\serion

Page 2/8

SERION ELISA classic Kit

SECTION 4: First aid measures

4.1. Description of first aid measures

General information:

In case of accident or unwellness, seek medical advice immediately (show directions for use or safety data sheet if possible). Remove victim out of the danger area. Remove contaminated, saturated clothing. If unconscious but breathing normally, place in recovery position and seek medical advice. Do not leave affected person unattended.

Following inhalation:

Provide fresh air.

After eye contact:

Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

Following ingestion:

Rinse mouth. Let 1 glass of water be drunken in little sips (dilution effect). Get medical advice/attention if you feel unwell.

4.2. Most important symptoms and effects, both acute and delayed

No known symptoms to date.

4.3. Indication of any immediate medical attention and special treatment needed

Treat symptomatically.

SECTION 5: Firefighting measures

5.1. Extinguishing media

No data available

5.2. Special hazards arising from the substance or mixture

Hazardous combustion products:

In case of fire: Gases/vapours, toxic

5.3. Advice for firefighters

Wear a self-contained breathing apparatus and chemical protective clothing.

5.4. Additional information

Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Personal precautions:

Avoid breathing dust/fume/gas/mist/vapours/spray. Remove persons to safety.

Protective equipment:

Wear protective gloves/protective clothing/eye protection/face protection.

6.1.2. For emergency responders

Personal protection equipment:

Personal protection equipment: see section 8

6.2. Environmental precautions

Do not allow to enter into surface water or drains.

6.3. Methods and material for containment and cleaning up

For containment:

Collect spillage. Measures to prevent aerosol and dust generation Wet clean or vacuum up solids.

For cleaning up:

Water (with cleaning agent)

6.4. Reference to other sections

Safe handling: see section 7 Personal protection equipment: see section 8 Disposal: see section 13

SAFETY INFORMATION

Voluntary safety information following the Safety Data Sheet format according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 26 Mar 2024

Print date: 26 Mar 2024

Version: 1.1

virion\serion

Page 3/8

SERION ELISA classic Kit

6.5. Additional information

Use appropriate container to avoid environmental contamination.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Protective measures

Advices on safe handling:

Wear personal protection equipment (refer to section 8).

Measures to prevent aerosol and dust generation:

Dust should be exhausted directly at the point of origin.

Advices on general occupational hygiene

When using do not eat, drink, smoke, sniff. Avoid contact with skin, eyes and clothes.

7.2. Conditions for safe storage, including any incompatibilities

Technical measures and storage conditions:

Keep container tightly closed in a cool, well-ventilated place.

7.3. Specific end use(s)

No data available

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

8.1.1. Occupational exposure limit values

Limit value type (country of origin)	Substance name	① Long-term occupational exposure limit value ② Short-term occupational exposure limit value ③ Instantaneous value ④ Monitoring and observation processes ⑤ Remark
WEL (GB)	sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1	① 0.1 mg/m ³ ② 0.3 mg/m ³ ⑤ (may be absorbed through the skin)
IOELV (EU)	sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1	① 0.1 mg/m ³ ② 0.3 mg/m ³ ⑤ (may be absorbed through the skin)
WEL (GB)	sodium hydroxide CAS No.: 1310-73-2 EC No.: 215-185-5	② 2 mg/m ³
IOELV (EU)	hydrogen chloride CAS No.: 7647-01-0 EC No.: 231-595-7	① 5 ppm (8 mg/m ³) ② 10 ppm (15 mg/m ³) ⑤ (Hydrogen chloride)
WEL (GB)	hydrogen chloride CAS No.: 7647-01-0 EC No.: 231-595-7	① 1 ppm (2 mg/m ³) ② 5 ppm (8 mg/m ³) ⑤ (gas and aerosol mists)

8.1.2. Biological limit values

No data available

SAFETY INFORMATION

Voluntary safety information following the Safety Data Sheet format according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 26 Mar 2024

Print date: 26 Mar 2024

Version: 1.1

virion\serion

Page 4/8

SERION ELISA classic Kit

8.1.3. DNEL-/PNEC-values

Substance name	DNEL value	① DNEL type ② Exposure route
sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1	0.0164	① DNEL worker ② Long-term - inhalation, systemic effects
sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1	0.0467	① DNEL worker ② Long-term - dermal, systemic effects
sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1	0.0167	① DNEL worker ② Long-term - oral, systemic effects

Substance name	PNEC Value	① PNEC type
Ethylendiamin-tetraessigsäure Dinatriumsalz CAS No.: 6381-92-6 EC No.: 205-358-3	2.2 mg/L	① PNEC aquatic, freshwater
Ethylendiamin-tetraessigsäure Dinatriumsalz CAS No.: 6381-92-6 EC No.: 205-358-3	0.22 mg/L	① PNEC aquatic, marine water
Ethylendiamin-tetraessigsäure Dinatriumsalz CAS No.: 6381-92-6 EC No.: 205-358-3	43 mg/L	① PNEC sewage treatment plant
sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1	0.00035 mg/L	① PNEC aquatic, freshwater
sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1	0.000015 mg/ L	① PNEC aquatic, marine water
sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1	0.0167 mg/kg	① PNEC sediment, freshwater
sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1	0.00072 mg/ kg	① PNEC sediment, marine water

8.2. Exposure controls

8.2.1. Appropriate engineering controls

No data available

8.2.2. Personal protection equipment

Eye/face protection:

Eye glasses with side protection EN 166

Skin protection:

Tested protective gloves must be worn EN ISO 374 Suitable material: Breakthrough time: min In the case of wanting to use the gloves again, clean them before taking off and air them well.

Respiratory protection:

Particle filter device (EN 143)

8.2.3. Environmental exposure controls

No data available

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Appearance

Physical state: solid

Colour: not determined

SAFETY INFORMATION

Voluntary safety information following the Safety Data Sheet format according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 26 Mar 2024

Print date: 26 Mar 2024

Version: 1.1

virion\serion

Page 5/8

SERION ELISA classic Kit

Odour: not determined

Safety relevant basis data

Parameter	Value	① Method ② Remark
pH	No data available	
Melting point	No data available	
Freezing point	No data available	
Initial boiling point and boiling range	No data available	
Evaporation rate	No data available	
Vapour pressure	No data available	
Density	No data available	
Bulk density	No data available	
Water solubility	No data available	

particle characteristics:

No data available

9.2. Other information

No data available

SECTION 10: Stability and reactivity

10.1. Reactivity

No data available

10.2. Chemical stability

No data available

10.3. Possibility of hazardous reactions

No data available

10.4. Conditions to avoid

No data available

10.5. Incompatible materials

No data available

10.6. Hazardous decomposition products

No data available

SECTION 11: Toxicological information

11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1

LD₅₀ oral: 27 mg/kg (Rat)

LD₅₀ dermal: 18 mg/kg (Rabbit)

LC₅₀ Acute inhalation toxicity (dust/mist): 5.4 mg/L 4 h (Rat)

Acute oral toxicity:

Based on available data, the classification criteria are not met.

Acute dermal toxicity:

Based on available data, the classification criteria are not met.

Acute inhalation toxicity:

Based on available data, the classification criteria are not met.

Skin corrosion/irritation:

Based on available data, the classification criteria are not met.

Serious eye damage/irritation:

Based on available data, the classification criteria are not met.

SAFETY INFORMATION

Voluntary safety information following the Safety Data Sheet format according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 26 Mar 2024

Print date: 26 Mar 2024

Version: 1.1

virion\serion

Page 6/8

SERION ELISA classic Kit

Respiratory or skin sensitisation:

Contains 2-methyl-2H-isothiazol-3-one. May produce an allergic reaction.

Germ cell mutagenicity:

Based on available data, the classification criteria are not met.

Carcinogenicity:

Based on available data, the classification criteria are not met.

Reproductive toxicity:

Based on available data, the classification criteria are not met.

STOT-single exposure:

Based on available data, the classification criteria are not met.

STOT-repeated exposure:

Based on available data, the classification criteria are not met.

Aspiration hazard:

Based on available data, the classification criteria are not met.

Additional information:

No data available

11.2. Information on other hazards

No data available

SECTION 12: Ecological information

12.1. Toxicity

sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1

EC₅₀: 0.348 mg/L 4 d (Algae/water plant)

EC₅₀: 0.348 mg/L 4 d

12.2. Persistence and degradability

sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1

Biodegradation: —

12.3. Bioaccumulative potential

No data available

12.4. Mobility in soil

No data available

12.5. Results of PBT and vPvB assessment

sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1

Results of PBT and vPvB assessment: —

hydrogen chloride CAS No.: 7647-01-0 EC No.: 231-595-7

Results of PBT and vPvB assessment: —

2-methyl-2H-isothiazol-3-one CAS No.: 2682-20-4 EC No.: 220-239-6

Results of PBT and vPvB assessment: —

12.6. Endocrine disrupting properties

No data available

12.7. Other adverse effects

No data available

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment options

Appropriate disposal / Product:

Consult the appropriate local waste disposal expert about waste disposal.

SAFETY INFORMATION

Voluntary safety information following the Safety Data Sheet format according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 26 Mar 2024

Print date: 26 Mar 2024

Version: 1.1

virion\serion

Page 7/8

SERION ELISA classic Kit

SECTION 14: Transport information

Land transport (ADR/RID)	Inland waterway craft (ADN)	Sea transport (IMDG)	Air transport (ICAO-TI / IATA-DGR)
14.1. UN number or ID number			
No dangerous good in sense of these transport regulations.	No dangerous good in sense of these transport regulations.	No dangerous good in sense of these transport regulations.	No dangerous good in sense of these transport regulations.
14.2. UN proper shipping name			
No dangerous good in sense of these transport regulations.	No dangerous good in sense of these transport regulations.	No dangerous good in sense of these transport regulations.	No dangerous good in sense of these transport regulations.
14.3. Transport hazard class(es)			
not relevant	not relevant	not relevant	not relevant
14.4. Packing group			
not relevant	not relevant	not relevant	not relevant
14.5. Environmental hazards			
not relevant	not relevant	not relevant	not relevant
14.6. Special precautions for user			
not relevant	not relevant	not relevant	not relevant

14.7. Maritime transport in bulk according to IMO instruments

No data available

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

No data available

15.2. Chemical Safety Assessment

No data available

SECTION 16: Other information

* 16.1. Indication of changes

1.4.	Emergency telephone number
16.1.	Indication of changes

16.2. Abbreviations and acronyms

ACGIH	American Conference of Governmental Industrial Hygienists
ADN	European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
CAS	Chemical Abstracts Service
CLP	Classification, Labelling and Packaging
DIN	German Institute for Standardization / German Industrial Standard
DNEL	derived no-effect level
EC ₅₀	Effective Concentration 50%
EN	European Standard
ES	Exposure scenario
ICAO	International Civil Aviation Organization
IMDG	International Maritime Dangerous Goods
IMO	International Maritime Organization
ISO	International Standards Organisation
LC ₅₀	Lethal (fatal) Concentration 50%
LD ₅₀	Lethal (fatal) Dose 50%
MAK	Maximum concentration in the workplace air (CH)

SAFETY INFORMATION

Voluntary safety information following the Safety Data Sheet format according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 26 Mar 2024

Print date: 26 Mar 2024

Version: 1.1

virion\serion

Page 8/8

SERION ELISA classic Kit

NFPA	National Fire Protection Association
NIOSH	National Institute for Occupational Safety & Health
OSHA	Occupational Safety & Health Administration
PBT	persistent and bioaccumulative and toxic
PNEC	Predicted No Effect Concentration
REACH	Registration, Evaluation and Authorization of Chemicals
RID	Dangerous goods regulations for transport by rail
SCL	Specific concentration limit
TRGS	Technische Regeln für Gefahrstoffe
UN	United Nations

16.3. Key literature references and sources for data

No data available

16.4. Classification for mixtures and used evaluation method according to regulation (EC) No 1272/2008 [CLP]

The mixture is classified as not hazardous according to regulation (EC) No 1272/2008 [CLP].

16.5. List of relevant hazard statements and/or precautionary statements from sections 2 to 15

No data available

16.6. Training advice

No data available

16.7. Additional information

No data available

* Data changed compared with the previous version.