Voluntary safety information following the Safety Data Sheet format according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 5 Apr 2024 Print date: 5 Apr 2024 Version: 1

virion\serion

Page 1/8

SERION ELISA Avidity Reagent

SECTION 1: Identification of the substance/mixture and of the company/ undertaking

1.1. Product identifier Trade name/designation:

SERION ELISA Avidity Reagent

Article No.:

B*AVID

1.2. Relevant identified uses of the substance or mixture and uses advised against Use of the substance/mixture:

The product is intended for professional use.

1.3. Details of the supplier of the safety data sheet

Supplier (manufacturer/importer/only representative/downstream user/distributor): Institut Virion\Serion GmbH

Productmanagement Friedrich-Bergius-Ring 19 97076 Würzburg Germany **Telephone:** +49 931 3045 0 **Telefax:** +49 931 3045 100 **E-mail:** product.safety@virion-serion.de **Website:** www.virion-serion.de

E-mail (competent person): product.safety@virion-serion.de

1.4. Emergency telephone number

Institut Virion\Serion GmbH, Friedrich-Bergius-Ring 19; 97076 Würzburg, Germany, +49 931 3045 0 (Only available during office hours.)

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No 1272/2008 [CLP]

The mixture is classified as not hazardous according to regulation (EC) No 1272/2008 [CLP].

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

According to EC directives or the corresponding national regulations the product does not have to be labelled.

Hazard components for labelling:

Citronensäure Monohydrat

Hazard statements: none

Supplemental hazard information: none

Precautionary statements: none

2.3. Other hazards

No data available

SAFETY INFORMATION Voluntary safety information following the Safety Data Sheet format according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 5 Apr 2024 Print date: 5 Apr 2024 Version: 1

virion\serion

Page 2/8

SERION ELISA Avidity Reagent

SECTION 3: Composition/information on ingredients

3.2. Mixtures

Hazardous ingredients / Hazardous impurities / Stabilisers:

Product identifiers	Substance name Classification according to Regulation (EC) No 1272/2008 [CLP]	Concentration
CAS No.: 5949-29-1 EC No.: 611-842-9 REACH No.: 01-2119457026-42-XXXX	Citronensäure Monohydrat Eye Irrit. 2 (H319), STOT SE 3 (H335) Warning	2 - ≤ 5 weight-%
CAS No.: 26628-22-8 EC No.: 247-852-1 Full text of H- and EUH-phras	sodium azide Substance with a community workplace exposure limit. ses: see section 16.	0 - ≤ 0.1 weight-%

SECTION 4: First aid measures

4.1. Description of first aid measures

General information:

In case of accident or unwellness, seek medical advice immediately (show directions for use or safety data sheet if possible). Remove victim out of the danger area. Remove contaminated, saturated clothing. If unconscious but breathing normally, place in recovery position and seek medical advice. Do not leave affected person unattended.

Following inhalation:

Provide fresh air.

After eye contact:

Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing. If eye irritation persists: Get medical advice/attention.

Following ingestion:

Rinse mouth. Let 1 glass of water be drunken in little sips (dilution effect). Get medical advice/attention if you feel unwell.

4.2. Most important symptoms and effects, both acute and delayed No known symptoms to date.

4.3. Indication of any immediate medical attention and special treatment needed Treat symptomatically.

SECTION 5: Firefighting measures

5.1. Extinguishing media

No data available

5.2. Special hazards arising from the substance or mixture

Hazardous combustion products:

In case of fire: Gases/vapours, toxic

5.3. Advice for firefighters

Wear a self-contained breathing apparatus and chemical protective clothing.

5.4. Additional information

Collect contaminated fire extinguishing water separately. Do not allow entering drains or surface water.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Personal precautions:

Avoid breathing dust/fume/gas/mist/vapours/spray. Remove persons to safety.

Protective equipment:

Wear protective gloves/protective clothing/eye protection/face protection.

Voluntary safety information following the Safety Data Sheet format according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 5 Apr 2024 Print date: 5 Apr 2024 Version: 1

virion**serion**

Page 3/8

SERION ELISA Avidity Reagent

6.1.2. For emergency responders

Personal protection equipment:

Personal protection equipment: see section 8

6.2. Environmental precautions

Do not allow to enter into surface water or drains.

6.3. Methods and material for containment and cleaning up

For containment:

Collect spillage. Measures to prevent aerosol and dust generation Wet clean or vacuum up solids.

For cleaning up:

Water (with cleaning agent)

6.4. Reference to other sections

Safe handling: see section 7 Personal protection equipment: see section 8 Disposal: see section 13

6.5. Additional information

Use appropriate container to avoid environmental contamination.

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Protective measures

Advices on safe handling:

Wear personal protection equipment (refer to section 8).

Measures to prevent aerosol and dust generation:

Dust should be exhausted directly at the point of origin.

Advices on general occupational hygiene

When using do not eat, drink, smoke, sniff. Avoid contact with skin, eyes and clothes.

7.2. Conditions for safe storage, including any incompatibilities

Technical measures and storage conditions:

Keep container tightly closed in a cool, well-ventilated place.

7.3. Specific end use(s)

No data available

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

8.1.1. Occupational exposure limit values

Limit value type (country of origin)	Substance name	 Long-term occupational exposure limit value Short-term occupational exposure limit value Instantaneous value Monitoring and observation processes Remark
WEL (GB)	sodium hydroxide CAS No.: 1310-73-2 EC No.: 215-185-5	② 2 mg/m ³
WEL (GB)	sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1	 0.1 mg/m³ 0.3 mg/m³ (may be absorbed through the skin)
IOELV (EU)	sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1	 0.1 mg/m³ 0.3 mg/m³ (may be absorbed through the skin)

8.1.2. Biological limit values

No data available

SAFETY INFORMATION Voluntary safety information following the Safety Data Sheet format according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 5 Apr 2024 Print date: 5 Apr 2024 Version: 1

virion\serion

Page 4/8

SERION ELISA Avidity Reagent

8.1.3. DNEL-/PNEC-values

Substance name	stance name DNEL value 🗓 DNEL type	
		② Exposure route
sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1	0.0164	 DNEL worker Long-term - inhalation, systemic effects
sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1	0.0467	 DNEL worker Long-term - dermal, systemic effects
sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1	0.0167	 DNEL worker Long-term - oral, systemic effects
Substance name	PNEC Value	① PNEC type
sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1	0.00035 mg/L	① PNEC aquatic, freshwater
sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1	0.000015 mg/ L	${f 1}$ PNEC aquatic, marine water
sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1	0.0167 mg/kg	${f 1}$ PNEC sediment, freshwater
sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1	0.00072 mg/ kg	 PNEC sediment, marine water

8.2. Exposure controls

8.2.1. Appropriate engineering controls

No data available

8.2.2. Personal protection equipment

Eye/face protection:

Eye glasses with side protection EN 166

Skin protection:

Tested protective gloves must be worn EN ISO 374 Suitable material: NBR (Nitrile rubber) Breakthrough time: 480 min In the case of wanting to use the gloves again, clean them before taking off and air them well. Breakthrough times and swelling properties of the material must be taken into consideration.

Respiratory protection:

Particle filter device (EN 143)

8.2.3. Environmental exposure controls

No data available

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Appearance

Physical state: solid Odour: not determined Colour: not determined

Odour: not determined Safety relevant basis data

Parameter	Value	 Method Remark
рН	No data available	
Melting point	No data available	
Freezing point	No data available	
Initial boiling point and boiling range	No data available	
Flash point	No data available	

Voluntary safety information following the Safety Data Sheet format according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 5 Apr 2024 Print date: 5 Apr 2024 Version: 1

virion\serion

Page 5/8

SERION ELISA Avidity Reagent

Parameter	Value	 Method 	
		② Remark	
Evaporation rate	No data available		
Auto-ignition temperature	No data available		
Upper/lower flammability or explosive limits	No data available		
Vapour pressure	No data available		
Vapour density	No data available		
Density	No data available		
Bulk density	not applicable		
Water solubility	No data available		
Dynamic viscosity	No data available		
Kinematic viscosity	No data available		

9.2. Other information

No data available

SECTION 10: Stability and reactivity

10.1. Reactivity

No data available

10.2. Chemical stability

No data available

10.3. Possibility of hazardous reactions No data available

10.4. Conditions to avoid

No data available

10.5. Incompatible materials

No data available

10.6. Hazardous decomposition products

No data available

SECTION 11: Toxicological information

11.1. Information on hazard classes as defined in Regulation (EC) No 1272/2008

Citronensäure Monohydrat CAS No.: 5949-29-1 EC No.: 611-842-9

LD₅₀ oral: 5,400 mg/kg (Mouse)

LD₅₀ dermal: >2,000 mg/kg (Rat)

sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1

LD₅₀ oral: 27 mg/kg (Rat)

LD₅₀ dermal: 18 mg/kg (Rabbit)

LC₅₀ Acute inhalation toxicity (dust/mist): 5.4 mg/L 4 h (Rat)

Acute oral toxicity:

Based on available data, the classification criteria are not met.

Acute dermal toxicity:

Based on available data, the classification criteria are not met.

Acute inhalation toxicity:

Based on available data, the classification criteria are not met.

Skin corrosion/irritation:

Based on available data, the classification criteria are not met.

Serious eye damage/irritation:

Based on available data, the classification criteria are not met.

Voluntary safety information following the Safety Data Sheet format according to Regulation (EC) No. 1907/2006 (REACH) **Revision date:** 5 Apr 2024 **Print date:** 5 Apr 2024 **Version:** 1

virion**serio**n

Page 6/8

SERION ELISA Avidity Reagent

Respiratory or skin sensitisation:

Based on available data, the classification criteria are not met.

Germ cell mutagenicity: Based on available data, the classification criteria are not met.

Carcinogenicity:

Based on available data, the classification criteria are not met.

Reproductive toxicity:

Based on available data, the classification criteria are not met.

STOT-single exposure:

Based on available data, the classification criteria are not met.

STOT-repeated exposure:

Based on available data, the classification criteria are not met.

Aspiration hazard:

Based on available data, the classification criteria are not met.

Additional information:

No data available

11.2. Information on other hazards No data available

SECTION 12: Ecological information

12.1. Toxicity

Citronensäure Monohydrat CAS No.: 5949-29-1 EC No.: 611-842-9

LC₅₀: 440 mg/L 2 d (fish)

sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1

EC₅₀: 0.348 mg/L 4 d (Algae/water plant)

EC₅₀: 0.348 mg/L 4 d

12.2. Persistence and degradability

sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1

Biodegradation: —

12.3. Bioaccumulative potential

No data available

12.4. Mobility in soil

No data available

12.5. Results of PBT and vPvB assessment

Citronensäure Monohydrat CAS No.: 5949-29-1 EC No.: 611-842-9

Results of PBT and vPvB assessment: -

sodium azide CAS No.: 26628-22-8 EC No.: 247-852-1

Results of PBT and vPvB assessment: -

12.6. Endocrine disrupting properties

No data available

12.7. Other adverse effects

No data available

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment options

Appropriate disposal / Product:

Consult the appropriate local waste disposal expert about waste disposal.

SAFETY INFORMATION Voluntary safety information following the Safety Data Sheet format according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 5 Apr 2024 Print date: 5 Apr 2024 Version: 1

virion**serion**

Page 7/8

SERION ELISA Avidity Reagent

SECTION 14: Transport information

Land transport (ADR/RID) Inland waterway craft (ADN)	Sea transport (IMDG)	Air transport (ICAO-TI / IATA-DGR)
14.1. UN number or	ID number		
No dangerous good in sense of these transport regulations.	No dangerous good in sense of these transport regulations.	No dangerous good in sense of these transport regulations.	No dangerous good in sense of these transport regulations.
14.2. UN proper shi	pping name	•	
No dangerous good in sense of these transport regulations.	No dangerous good in sense of these transport regulations.	No dangerous good in sense of these transport regulations.	No dangerous good in sense of these transport regulations.
14.3. Transport haz	ard class(es)	·	
not relevant	not relevant	not relevant	not relevant
14.4. Packing group)		
not relevant	not relevant	not relevant	not relevant
14.5. Environmenta	l hazards		
not relevant	not relevant	not relevant	not relevant
14.6. Special preca	utions for user		
not relevant	not relevant	not relevant	not relevant

14.7. Maritime transport in bulk according to IMO instruments No data available

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

No data available

15.2. Chemical Safety Assessment

No data available

SECTION 16: Other information

16.1. Indication of changes

No data available

16.2. Abbreviations and acronyms

- ACGIH American Conference of Governmental Industrial Hygienists
- ADN European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways
- ADR European Agreement concerning the International Carriage of Dangerous Goods by Road CAS Chemical Abstracts Service
- CLP Classification, Labelling and Packaging
- DIN German Institute for Standardization / German Industrial Standard
- DNEL derived no-effect level
- EC₅₀ Effective Concentration 50%
- EN European Standard
- ES Exposure scenario
- ICAO International Civil Aviation Organization
- IMDG International Maritime Dangerous Goods
- IMO International Maritime Organization
- ISO International Standards Organisation
- LC₅₀ Lethal (fatal) Concentration 50%
- LD₅₀ Lethal (fatal) Dose 50%
- MAK Maximum concentration in the workplace air (CH)
- NFPA National Fire Protection Association

Voluntary safety information following the Safety Data Sheet format according to Regulation (EC) No. 1907/2006 (REACH)

Revision date: 5 Apr 2024 Print date: 5 Apr 2024 Version: 1

virion**serion**

Page 8/8

SERION ELISA Avidity Reagent

- NIOSH National Institute for Occupational Safety & Health
- OSHA Occupational Safety & Health Administration
- PBT persistent and bioaccumulative and toxic
- PNEC Predicted No Effect Concentration
- REACH Registration, Evaluation and Authorization of Chemicals
- RID Dangerous goods regulations for transport by rail
- TRGS Technische Regeln für Gefahrstoffe
- UN United Nations

16.3. Key literature references and sources for data

No data available

16.4. Classification for mixtures and used evaluation method according to regulation (EC) No 1272/2008 [CLP]

The mixture is classified as not hazardous according to regulation (EC) No 1272/2008 [CLP].

16.5. List of relevant hazard statements and/or precautionary statements from sections 2 to 15

Hazard statements	
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.

16.6. Training advice

No data available

16.7. Additional information

No data available

en / GB GeSi.de