

# Analysenzertifikat

## Certificate of Analysis

### SERION ELISA *classic*

### Zusatzreagenz / Supplementary Reagent

### SERION RF Absorbent

Best.-Nr. / Order No.	Z200
Chargen-Nr. / Lot No.	CQ0106
IFU-Version / IFU version	12-25/06-2
Herstelldatum / Manufacturing date	2025-07
Verfallsdatum / Expiry Date	2027-07-31
Anwendungsgebiet / Field of Application	Entfernung von Rheumafaktoren durch Absorption von 15 mg/ml IgG Removal of rheumatoid factor by absorption of 15 mg/ml IgG
Material / Material	Ziegenserum / Goat serum
Konservierungsmittel / Preservation	0,09% NaN <sub>3</sub>
Lagerung / Storage	2°C - 8°C
Gebrauchsverdünnung / Working Dilution	gebrauchsfertig / ready to use



Hersteller / Manufacturer:  
 Institut Virion\Serion GmbH  
 Friedrich-Bergius-Ring 19  
 D-97076 Würzburg  
 Germany  
 Tel. +49 931 3045 0  
 Fax +49 931 3045 100  
 E-Mail: [info@virion-serion.de](mailto:info@virion-serion.de)  
 Internet: [www.virion-serion.de](http://www.virion-serion.de)