

Qualitätskontrollzertifikat

Quality Control Certificate

SERION ELISA *classic* Coxiella burnetii Phase 1 IgA

Best.-Nr. / Order No.	ESR1311A
Chargen-Nr. / Lot No.	EQ0216
Verfallsdatum / Expiry Date	2027-10-31
IFU-Version / IFU-Version	NEW! 131-17

Kontrollseren / Control sera	Chargen-Nr. / Lot No	Gemessener Wert / Measured value OD 405 nm	Gültigkeitsbereich / Validity range
Positivkontrolle / Positive control	ECQ0473	1,31	> 0,66
Cut off-Serum / Cut off serum	ECQ0474	0,41	0,21 – 0,82

Negativkontrolle / Negative Control Chargen-Nr. / Lot No.	Coxiella burnetii Phase 1 IgA ECQ0472
Teststreifen / Strips Chargen-Nr. / Lot No.	Coxiella burnetii Phase 1 IgA ECQ0468
Konjugat / Conjugate Chargen-Nr. / Lot No.	IgA ++ KJQ093++

31.10.2025

Prüfdatum / Date of control

Hersteller / Manufacturer:
 Institut Virion\Serion GmbH
 Friedrich-Bergius-Ring 19
 D-97076 Würzburg
 Germany
 Tel. +49 931 3045 0
 Fax +49 931 3045 100
 E-Mail: info@virion-serion.de
 Internet: www.virion-serion.de