

Qualitätskontrollzertifikat

Quality Control Certificate

SERION ELISA *classic* Coxiella burnetii Phase 2 IgM

Best.-Nr. / Order No.	ESR1312M
Chargen-Nr. / Lot No.	ER0106
Verfallsdatum / Expiry Date	2028-01-31
IFU-Version / IFU-Version	NEW! 131-17

Kontrollseren / Control sera	Chargen-Nr. / Lot No	Gemessener Wert / Measured value OD 405 nm	Gültigkeitsbereich / Validity range
Positivkontrolle / Positive control	ECR0143	0,98	> 0,49
Cut off-Serum / Cut off serum	ECR0186	0,37	0,19 - 0,74

Negativkontrolle / Negative Control	Coxiella burnetii Phase 2 IgM
Chargen-Nr. / Lot No.	ECR0142
Teststreifen / Strips	Coxiella burnetii Phase 2 IgM
Chargen-Nr. / Lot No.	ECQ0146
Konjugat / Conjugate	IgM +
Chargen-Nr. / Lot No.	KJQ090

27.04.2026

Prüfdatum / Date of control

Hersteller / Manufacturer:
 Institut Virion\Serion GmbH
 Friedrich-Bergius-Ring 19
 D-97076 Würzburg
 Germany
 Tel. +49 931 3045 0
 Fax +49 931 3045 100
 E-Mail: info@virion-serion.de
 Internet: www.virion-serion.de