

Analysenzertifikat**Certificate of Analysis****SERION *immuno* CONTROL (ToRCH IgG)****Toxoplasma gondii IgG****Rubella Virus IgG****Cytomegalovirus IgG****Herpes Simplex Virus 1+2 IgG**

Artikelnr. / Article No.	IC1701
Chargen-Nr. / Lot No.	SCI.DL
Herstelldatum / Manufacturing Date	2018-03
Verfallsdatum / Expiry Date	2020-03
Versions-Nr. / Version No.	2018 / 03 - 1

zur Verwendung mit dem entsprechenden SERION ELISA *classic* Test:
for use with the corresponding SERION ELISA *classic* assay:

Rekonstitution / Reconstitution	1 ml <i>Aqua dest.</i>
Zielwert / Target value	positiv / positive



Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de