

# Qualitätskontrollzertifikat

## Quality Control Certificate

### SERION ELISA *classic* ADENOVIRUS IgM

Best.-Nr. / Order No.	ESR128M
Chargen-Nr. / Lot No.	EK0068
Verfallsdatum / Expiry Date	2021-04
IFU-Version / IFU-Version	<b>NEW!</b> 128-15

Kontrollseren / Control sera	Chargen-Nr. / Lot No	Gemessener Wert / Measured value OD 405 nm	Gültigkeitsbereich / Validity range
Positivkontrolle / Positive control	ECK0079	1,01	> 0,51
Cut off-Serum / Cut off serum	ECK0080	0,26	0,13 – 0,52

Negativkontrolle / Negative Control	ADENOVIRUS IgM
Chargen-Nr. / Lot No.	ECK0078
Teststreifen / Strips	ADENOVIRUS IgM
Chargen-Nr. / Lot No.	SDK.DY
Konjugat / Conjugate	IgM ++
Chargen-Nr. / Lot No.	KJK002

15.07.2019

Prüfdatum / Date of control

Hersteller / Manufacturer:  
 Institut Virion\Serion GmbH  
 Friedrich-Bergius-Ring 19  
 D-97076 Würzburg  
 Germany  
 Tel. +49 931 3045 0  
 Fax +49 931 3045 100  
 E-Mail: info@virion-serion.de  
 Internet: www.virion-serion.de