

Analysenzertifikat

Certificate of Analysis

SERION ELISA control Aviditätstest / Avidity Test CYTOMEGALOVIRUS

Best.-Nr. / Order No.	BR109AVID
Chargen-Nr. / Lot No.	SDI.CQ
Herstelldatum / Manufacturing Date	2018-04
Verfallsdatum / Expiry Date	2028-04
Zielwert (SERION Aviditätsindex) / Target Value (SERION Avidity Index)	Zielwert 29,6%
Zielbereich (SERION Aviditätsindex) / Target Range (SERION Avidity Index)	$\geq 14,8 \%$ und / and $\leq 44,4 \%$
Rekonstitution / Reconstitution	0,1 ml Aqua dest.
Serumverdünnung / Serum Dilution	1 + 100
Versions-Nr. / Version No.	2018 / 04-1



Hersteller / Manufacturer:

Institut Virion\Serion GmbH

Friedrich-Bergius-Ring 19

D-97076 Würzburg

Germany

Tel. +49 931 3045 0

Fax +49 931 3045 100

E-Mail: dialog@virion-serion.de

Internet: www.virion-serion.de