



Analysenzertifikat Certificate of Analysis

virion\sersion

SERION ELISA *classic* Aviditätsreagenz / Avidity Reagent CYTOMEGALOVIRUS

Best.-Nr. / Order No.	B109AVID
Chargen-Nr. / Lot No.	SHD.BR
Verfallsdatum / Expiry Date	2023-08
Verwendung mit / For use with SERION ELISA <i>classic</i> Cytomegalovirus IgG	SGD.FE
Verwendung mit Standardserum / For use with standard serum	SGD.FD
Standard Zielwert (SERION Aviditätsindex) / Standard Target Value (SERION Avidity Index)	75 %
Standard Zielbereich (SERION Aviditätsindex) / Standard Target Range (SERION Avidity Index)	≥ 60 % und / and ≤ 90 %
X Faktor / X factor	0,428
Y Faktor / Y factor	1,534
Rekonstitution / Reconstitution	1,5 ml Aqua dest.
Versions-Nr. / Version No.	2013 / 12-1

Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de