

Analysenzertifikat

Certificate of Analysis

SERION ELISA *control* Aviditätstest / Avidity Test RÖTELN / RUBELLA IGG

Best.-Nr. / Order No.	BR129AVID
Chargen-Nr. / Lot No.	SEI.CE
Herstelldatum / Manufacturing Date	2018-05
Verfallsdatum / Expiry Date	2028-05
Zielwert (SERION Aviditätsindex) / Target Value (SERION Avidity Index)	6,8 %
Zielbereich (SERION Aviditätsindex) / Target Range (SERION Avidity Index)	$\geq 3,4$ % und / and $\leq 10,2$ %
Rekonstitution / Reconstitution	0,1 ml Aqua dest.
Serumverdünnung / Serum Dilution	1 + 100
Versions-Nr. / Version No.	2018 / 05-1



Hersteller / Manufacturer:

Institut Virion\Serion GmbH

Friedrich-Bergius-Ring 19

D-97076 Würzburg

Germany

Tel. +49 931 3045 0

Fax +49 931 3045 100

E-Mail: dialog@virion-serion.de

Internet: www.virion-serion.de