

Analysenzertifikat

Certificate of Analysis

SERION ELISA *control* Aviditätstest / Avidity Test RÖTELN / RUBELLA IgG

Best.-Nr. / Order No.	BR129AVID
Chargen-Nr. / Lot No.	SEI.CE
Herstelldatum / Manufacturing Date	2018-05
Verfallsdatum / Expiry Date	2028-05
Zielwert (SERION Aviditätsindex) / Target Value (SERION Avidity Index)	9,8 %
Zielbereich (SERION Aviditätsindex) / Target Range (SERION Avidity Index)	$\geq 4,9 \%$ und / and $\leq 19,6 \%$
Rekonstitution / Reconstitution	0,1 ml Aqua dest.
Serumverdünnung / Serum Dilution	1 + 100
Versions-Nr. / Version No.	2018 / 07-1
gültig ab SERION ELISA classic Kitcharge / valid from SERION ELISA classic lot	SEI.BU



Hersteller / Manufacturer:
 Institut Virion\Serion GmbH
 Friedrich-Bergius-Ring 19
 D-97076 Würzburg
 Germany
 Tel. +49 931 3045 0
 Fax +49 931 3045 100
 E-Mail: dialog@virion-serion.de
 Internet: www.virion-serion.de