


ut Virion\Serion GmbH · Postfach 92 53 · 97092 Würzburg

Analysenzertifikat Certificate of Analysis



virion\serion^g

Tel. +49 (0) 931 3045 0
Fax +49 (0) 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de

SERION KBR Positivkontrolle / CFT Positive Control PARAINFLUENZA VIRUS 2

Best.-Nr. / Order No.	3117
Chargen-Nr. / Lot No.	SHC.CN
Verfallsdatum / Expiry Date	2022-08
Material / Material	Kaninchenserum/ rabbit serum
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	0,1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:10
KBR Titer / CFT Titre	1:80

serion



multianalyt™

Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de