



Analysenzertifikat Certificate of Analysis

virion\serion

SERION KBR Positivkontrolle / CFT Positive Control YERSINIA PSEUDOTUBERCOSIS

| | |
|---|-------------------------------|
| Best.-Nr. / Order No. | 3201 |
| Chargen-Nr. / Lot No. | SBE.AD |
| Verfallsdatum / Expiry Date | 2024-02 |
| Material / Material | Kaninchenserum / rabbit serum |
| Lagerung / Storage | 2 °C - 8 °C |
| Rekonstitution / Reconstitution | 0,1 ml Aqua dest. |
| Gebrauchsverdünnung / Working Dilution | 1:10 |
| KBR Titer / CFT Titre | 1:40 |

serion



multianalyt™

Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de