



Institut Virion\Serion GmbH · Postfach 92 53 · 97092 Würzburg

Analysenzertifikat Certificate of Analysis

virion\serion
 Tel. +49 (0) 931 3045 0
 Fax +49 (0) 931 3045 100
 E-Mail: dialog@virion-serion.de
 Internet: www.virion-serion.de

SERION KBR Positivkontrolle / CFT Positive Control ERREGERNAME

| | |
|---|---|
| Best.-Nr. / Order No. | 3224 |
| Chargen-Nr. / Lot No. | SGC.AZ |
| Verfallsdatum / Expiry Date | 2022-07 |
| Material / Material | Meerschweinchenserum/guinea pig serum verdünnt/diluted |
| Lagerung / Storage | 2 °C - 8 °C |
| Rekonstitution / Reconstitution | 0,1 ml Aqua dest. |
| Gebrauchsverdünnung / Working Dilution | 1:10 |
| KBR Titer / CFT Titre | 1:40 |

serion



multianalyt™

Hersteller / Manufacturer:
 Institut Virion\Serion GmbH
 Friedrich-Bergius-Ring 19
 D-97076 Würzburg
 Germany
 Tel. +49 931 3045 0
 Fax +49 931 3045 100
 E-Mail: dialog@virion-serion.de
 Internet: www.virion-serion.de