

# Analysenzertifikat

## Certificate of Analysis

### SERION KBR Antigen / CFT Antigen ADENOVIRUS

Best.-Nr. / Order No.	1121
Chargen-Nr. / Lot No.	SKH.DH
Verfallsdatum / Expiry Date	2027-10
Antigen / Antigen	Adenoviren / Adenoviruses 2, 3, 4, 7 (Adenoid 6, GB, RI-67, Gomen)
Kultursystem / Culture System	Verozellen / Vero cells
Präparation / Preparation	Mischung aus infizierten Zellen und Zellkulturmedium (konzentriert und gereinigt)  Mixture of infected cells and cell culture medium (concentrated and purified)
Inaktivierung / Inactivation	Hitzebehandlung / heat treatment
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:32



Hersteller / Manufacturer:  
 Institut Virion\Serion GmbH  
 Friedrich-Bergius-Ring 19  
 D-97076 Würzburg  
 Germany  
 Tel. +49 931 3045 0  
 Fax +49 931 3045 100  
 E-Mail: [dialog@virion-serion.de](mailto:dialog@virion-serion.de)  
 Internet: [www.virion-serion.de](http://www.virion-serion.de)