



Analysenzertifikat Certificate of Analysis

virion\serion

SERION KBR Antigen / CFT Antigen CAMPYLOBACTER FETUS SSP. INTESTINALIS

| | |
|---|--|
| Best.-Nr. / Order No. | 1207 |
| Chargen-Nr. / Lot No. | SHD.AV |
| Verfallsdatum / Expiry Date | 2023-08 |
| Antigen / Antigen | Campylobacter fetus ssp. intestinalis |
| Kultursystem / Culture System | Plattenkultur / cultivation on agar plates |
| Präparation / Preparation | beschallter Ganzzellextrakt sonicated cell lysate |
| Inaktivierung / Inactivation | Zerstörung der Zellen durch Ultraschall disruption of cells by sonication |
| Konservierung / Preservation | 0,01 % Thiomersal / 0.01 % thiomersal |
| Lagerung / Storage | 2 °C - 8 °C |
| Rekonstitution / Reconstitution | 1 ml Aqua dest. |
| Gebrauchsverdünnung / Working Dilution | 1:20 |

serion



multianalyt™

Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de