



Analysenzertifikat Certificate of Analysis

virion\serion

SERION KBR Antigen / CFT Antigen CAMPYLOBACTER FETUS SSP. INTESTINALIS

Best.-Nr. / Order No.	1207
Chargen-Nr. / Lot No.	SBH.CO
Verfallsdatum / Expiry Date	2027-02
Antigen / Antigen	Campylobacter fetus ssp. intestinalis
Kultursystem / Culture System	Plattenkultur / cultivation on agar plates
Präparation / Preparation	beschallter Ganzzellextrakt sonicated cell lysate
Inaktivierung / Inactivation	Zerstörung der Zellen durch Ultraschall disruption of cells by sonication
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:20



Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de