



Analysenzertifikat Certificate of Analysis

virion\sersion

SERION KBR Antigen / CFT Antigen CHLAMYDIA

Best.-Nr. / Order No.	1122
Chargen-Nr. / Lot No.	SAD.BQ
Verfallsdatum / Expiry Date	2023-01
Antigen / Antigen	Chlamydia gruppenspezifisches Antigen / Chlamydia group specific antigen
Kultursystem / Culture System	Verozellen / Vero cells
Präparation / Preparation	Gereinigte Elementar- und Retikularkörperchen von <i>Chlamydia abortus</i> und <i>trachomatis</i> / Purified elementary and reticular bodies of <i>Chlamydia abortus</i> and <i>trachomatis</i>
Inaktivierung / Inactivation	Hitzeinaktivierung / Heat treatment
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:6

serion



multianalyt™

Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de

FB-IV-047 A