



Analysenzertifikat Certificate of Analysis

virion\sersion

SERION KBR Antigen / CFT Antigen CHLAMYDIA

Best.-Nr. / Order No.	1122
Chargen-Nr. / Lot No.	SCH.AA
Verfallsdatum / Expiry Date	2027-03
Antigen / Antigen	Chlamydia gruppenspezifisches Antigen / Chlamydia group specific antigen
Kultursystem / Culture System	Verozellen / Vero cells
Präparation / Preparation	Gereinigte Elementar-und Retikularkörperchen von Chlamydia abortus und trachomatis / Purified elementary and reticular bodies of Chlamydia abortus and trachomatis
Inaktivierung / Inactivation	Hitzeinaktivierung / Heat treatment
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:6



Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de

FB-IV-047 B