



Analysenzertifikat Certificate of Analysis

virion\sersion

SERION KBR Antigen / CFT Antigen COXIELLA BURNETII PHASE 1

Best.-Nr. / Order No.	1227
Chargen-Nr. / Lot No.	SFH.AC
Verfallsdatum / Expiry Date	2027-06
Antigen / Antigen	Coxiella burnetii (Nine Mile)
Kultursystem / Culture System	bebrütete Hühnereier incubated hen's eggs
Präparation / Preparation	Coxiella burnetii Zellbestandteile (konzentriert und gereinigt) Coxiella burnetii cell components (concentrated and purified)
Inaktivierung / Inactivation	Behandlung mit 0,5 % Formalin/ treatment with 0.5 % formalin
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:8



Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de
FB-IV-047 B