



# Analysenzertifikat Certificate of Analysis

virion\serion

## SERION KBR Antigen / CFT Antigen COXIELLA BURNETII PHASE 2

Best.-Nr. / Order No.	1123
Chargen-Nr. / Lot No.	SIC.BB
Verfallsdatum / Expiry Date	2022-09
Antigen / Antigen	Coxiella burnetii (Nine Mile)
Kultursystem / Culture System	bebrütete Hühnereier incubated hen's eggs
Präparation / Preparation	Coxiella burnetii Zellbestandteile (konzentriert und gereinigt) Coxiella burnetii cell components (concentrated and purified)
Inaktivierung / Inactivation	Behandlung mit 0,5 % Formalin/ treatment with 0.5 % formalin
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:8

serion



multianalyt™

Hersteller / Manufacturer:  
Institut Virion\Serion GmbH  
Friedrich-Bergius-Ring 19  
D-97076 Würzburg  
Germany  
Tel. +49 931 3045 0  
Fax +49 931 3045 100  
E-Mail: [dialog@virion-serion.de](mailto:dialog@virion-serion.de)  
Internet: [www.virion-serion.de](http://www.virion-serion.de)