



# Analysenzertifikat Certificate of Analysis

virion\sersion

## SERION KBR Antigen / CFT Antigen COXSACKIEVIRUS A9

Best.-Nr. / Order No.	9060
Chargen-Nr. / Lot No.	SFD.BQ
Verfallsdatum / Expiry Date	2023-06
Antigen / Antigen	Coxsackievirus A9 Stamm / Strain P.B. (Bozek)
Kultursystem / Culture System	Verozellen / Vero cells
Präparation / Preparation	infizierte Zellen und Zellkulturmedium (konzentriert und gereinigt) infected cells and culture fluid (concentrated and purified)
Inaktivierung / Inactivation	beta-Propiolacton / beta-propiolactone
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:4

serion



multianalyt™

Hersteller / Manufacturer:  
Institut Virion\Serion GmbH  
Friedrich-Bergius-Ring 19  
D-97076 Würzburg  
Germany  
Tel. +49 931 3045 0  
Fax +49 931 3045 100  
E-Mail: [dialog@virion-serion.de](mailto:dialog@virion-serion.de)  
Internet: [www.virion-serion.de](http://www.virion-serion.de)