



Analysenzertifikat Certificate of Analysis

virion\sersion

SERION KBR Antigen / CFT Antigen COXSACKIE POOL (A9, B1-B6)

Best.-Nr. / Order No.	1178
Chargen-Nr. / Lot No.	SBF.FD
Verfallsdatum / Expiry Date	2025-02
Antigen / Antigen	Coxsackie Virus A9 und/and Coxsackie Virus B1 - B6
Kultursystem / Culture System	A9 - LLCMK2-Zellen/-cells B1 - B6 – Vero-Zellen/-cells
Präparation / Preparation	infizierte Zellen und Zellkulturmedium (konzentriert und gereinigt) infected cells and cell culture medium (concentrated and purified)
Inaktivierung / Inactivation	beta-Propiolacton und UV Behandlung / beta-propiolactone and UV treatment
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:4



Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de
FB-IV-047 B