



Analysenzertifikat Certificate of Analysis

virion\sersion

SERION KBR Antigen / CFT Antigen COXSACKIEVIRUS B6

Best.-Nr. / Order No.	1177
Chargen-Nr. / Lot No.	SGH.DT
Verfallsdatum / Expiry Date	2027-07
Antigen / Antigen	Coxsackievirus B6 (Schmitt)
Kultursystem / Culture System	Verozellen / vero cells
Präparation / Preparation	Infizierte Zellen und Zellkulturmedium (konzentriert und gereinigt) Infected cells and cell culture medium (concentrated and purified)
Inaktivierung / Inactivation	beta-Propiolacton / beta-propiolactone
Konservierung / Preservation	0,01 % Thiomersal / 0,01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:4



Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de
FB-IV-047 B