

# Analysenzertifikat

## Certificate of Analysis

### SERION KBR Antigen / CFT Antigen COXSACKIE POOL (A9, B1-B6)

|  |   |
|--|---|
| Best.-Nr. / Order No.                  | 1178  |
| Chargen-Nr. / Lot No.                  | SHH.CL  |
| Verfallsdatum / Expiry Date            | 2027-08   |
| Antigen / Antigen                      | Coxsackie Virus A9 und/and<br>Coxsackie Virus B1 - B6   |
| Kultursystem / Culture System          | A9 - LLCMK2-Zellen/-cells<br>B1 - B6 - Vero-Zellen/-cells   |
| Präparation / Preparation              | infizierte Zellen und Zellkulturmedium<br>(konzentriert und gereinigt)<br>infected cells and cell culture medium<br>(concentrated and purified) |
| Inaktivierung / Inactivation           | beta-Propiolacton und UV Behandlung /<br>beta-propiolactone and UV treatment  |
| Konservierung / Preservation           | 0,01 % Thiomersal / 0.01 % thiomersal   |
| Lagerung / Storage                     | 2 °C - 8 °C   |
| Rekonstitution / Reconstitution        | 1 ml Aqua dest.   |
| Gebrauchsverdünnung / Working Dilution | 1:4   |



Hersteller / Manufacturer:  
 Institut Virion\Serion GmbH  
 Friedrich-Bergius-Ring 19  
 D-97076 Würzburg  
 Germany  
 Tel. +49 931 3045 0  
 Fax +49 931 3045 100  
 E-Mail: [dialog@virion-serion.de](mailto:dialog@virion-serion.de)  
 Internet: [www.virion-serion.de](http://www.virion-serion.de)