



Analysenzertifikat Certificate of Analysis

virion\sersion

SERION KBR Antigen / CFT Antigen EPSTEIN-BARR VIRUS

Best.-Nr. / Order No.	1132
Chargen-Nr. / Lot No.	SCC.DP
Verfallsdatum / Expiry Date	2022-03
Antigen / Antigen	Epstein Barr Virus (EBV)
Kultursystem / Culture System	P3HR-1 Zelllinie / P3HR-1 cell line
Präparation / Preparation	Partiell gereinigte aufgeschlossene und alkalisch extrahierte Zellen / Partially purified disrupted and alkaline extracted cells
Inaktivierung / Inactivation	Inaktivierung durch alkalische Extraktion / Inactivation by alkaline extraction
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:4

serion



multianalyt™

Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de
FB-IV-047 A