



Analysenzertifikat Certificate of Analysis

virion\serion

SERION KBR Antigen / CFT Antigen ECHOVIRUS POOL

Best.-Nr. / Order No.	1180
Chargen-Nr. / Lot No.	SMD.CV
Verfallsdatum / Expiry Date	2023-12
Antigen / Antigen	Echoviren Typ / Echoviruses type 4, 6, 9, 14, 24, 30 (Stämme / strains: Pesascek, D'Amori, Hill, Tow, DeCamp, Bastianni)
Kultursystem / Culture System	Verozellen / Vero cells
Präparation / Preparation	infizierte Zellen und Zellkulturmedium (konzentriert und gereinigt) / infected cells and cell culture supernatant (concentrated and purified)
Inaktivierung / Inactivation	beta-Propiolacton und UV Behandlung / beta-propiolactone and UV treatment
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:6

serion



multianalyt™

Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de