



# Analysenzertifikat Certificate of Analysis

virion\sersion

## SERION KBR Antigen / CFT Antigen ECHOVIRUS POOL

Best.-Nr. / Order No.	1180
Chargen-Nr. / Lot No.	SIG.DH
Verfallsdatum / Expiry Date	2026-09
Antigen / Antigen	Echoviren Typ / Echoviruses type 4, 6, 9, 14, 24, 30 (Stämme / strains: Pesascek, D'Amori, Hill, Tow, DeCamp, Bastianni)
Kultursystem / Culture System	Verozellen / Vero cells
Präparation / Preparation	infizierte Zellen und Zellkulturmedium (konzentriert und gereinigt) / infected cells and cell culture supernatant (concentrated and purified)
Inaktivierung / Inactivation	beta-Propiolacton und UV Behandlung / beta-propiolactone and UV treatment
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:6



Hersteller / Manufacturer:  
Institut Virion\Serion GmbH  
Friedrich-Bergius-Ring 19  
D-97076 Würzburg  
Germany  
Tel. +49 931 3045 0  
Fax +49 931 3045 100  
E-Mail: [dialog@virion-serion.de](mailto:dialog@virion-serion.de)  
Internet: [www.virion-serion.de](http://www.virion-serion.de)

FB-IV-047 B