



Analysenzertifikat Certificate of Analysis

virion\serion

SERION KBR Antigen / CFT Antigen FSME VIRUS / TBE VIRUS

| | |
|---|--|
| Best.-Nr. / Order No. | 1192 |
| Chargen-Nr. / Lot No. | SCD.DF |
| Verfallsdatum / Expiry Date | 2023-03 |
| Antigen / Antigen | Frühsommer-Meningoenzephalitis- (FSME) Virus - Stamm Moskau B-4 / Tick borne encephalitis (TBE) virus - Strain Moscow B-4 |
| Kultursystem / Culture System | Verozellen / Vero cells |
| Präparation / Preparation | Partiell gereinigtes Virus aus Zellkultur- Überstand / Virus partially purified from cell culture supernatant |
| Inaktivierung / Inactivation | Behandlung mit β -Propiolacton / Treatment with β -propiolactone |
| Konservierung / Preservation | 0,01 % Thiomersal / 0.01 % thiomersal |
| Lagerung / Storage | 2 °C - 8 °C |
| Rekonstitution / Reconstitution | 1 ml Aqua dest. |
| Gebrauchsverdünnung / Working Dilution | 1:8 |

serion



multianalyt™