

# Analysenzertifikat

## Certificate of Analysis

**SERION KBR Antigen / CFT Antigen**  
**FSME VIRUS / TBE VIRUS**

Best.-Nr. / Order No.	1192
Chargen-Nr. / Lot No.	SLH.BL
Verfallsdatum / Expiry Date	2027-11
Antigen / Antigen	Frühsommer-Meningoenzephalitis-(FSME) Virus - Stamm Moskau B-4 / Tick borne encephalitis (TBE) virus - Strain Moscow B-4
Kultursystem / Culture System	Verozellen / Vero cells
Präparation / Preparation	Partiell gereinigtes Virus aus Zellkultur- Überstand / Virus partially purified from cell culture supernatant
Inaktivierung / Inactivation	Behandlung mit $\beta$ -Propiolacton / Treatment with $\beta$ -propiolactone
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:8



Hersteller / Manufacturer:  
 Institut Virion\Serion GmbH  
 Friedrich-Bergius-Ring 19  
 D-97076 Würzburg  
 Germany  
 Tel. +49 931 3045 0  
 Fax +49 931 3045 100  
 E-Mail: [dialog@virion-serion.de](mailto:dialog@virion-serion.de)  
 Internet: [www.virion-serion.de](http://www.virion-serion.de)