



Analysenzertifikat Certificate of Analysis

virion\serion

SERION KBR Antigen / CFT Antigen HERPES SIMPLEX VIRUS 1/2

Best.-Nr. / Order No.	1154
Chargen-Nr. / Lot No.	SHC.FD
Verfallsdatum / Expiry Date	2022-08
Antigen / Antigen	HSV1 (Stamm/strain F) HSV2 (Stamm/strain G)
Kultursystem / Culture System	Verozellen / Vero cell culture
Präparation / Preparation	Mischung aus infizierten Zellen und Zellkulturmedium (konzentriert und gereinigt) Mixture of infected cells and culture medium (concentrated and purified)
Inaktivierung / Inactivation	Hitzeinaktivierung und Tween-Ether Behandlung / heat inactivation and tween-ether treatment
Konservierung / Preservation	0,01 % Thiomersal 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:16

serion



multianalyt™