



Analysenzertifikat Certificate of Analysis

virion\serion

SERION KBR Antigen / CFT Antigen INFLUENZA B VIRUS

Best.-Nr. / Order No.	1113
Chargen-Nr. / Lot No.	SFE.AI
Verfallsdatum / Expiry Date	2024-06
Antigen / Antigen	Influenza B Virus (Stamm Brisbane) / Influenza B virus (strain Brisbane)
Kultursystem / Culture System	bebrütete Hühnereier / incubated hen's eggs
Präparation / Preparation	konzentrierte Virus Core-Partikel / concentrated virus core particles
Inaktivierung / Inactivation	Behandlung mit Detergens / treatment with detergent
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:16

serion



Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de