



# Analysenzertifikat Certificate of Analysis

virion\serion

## SERION KBR Antigen / CFT Antigen INFLUENZA B VIRUS

Best.-Nr. / Order No.	1113
Chargen-Nr. / Lot No.	SBF.CY
Verfallsdatum / Expiry Date	2025-02
Antigen / Antigen	Influenza B Virus (Stamm Brisbane)/ Influenza B Virus (strain Brisbane)
Kultursystem / Culture System	bebrütete Hühnereier / incubated hen's eggs
Präparation / Preparation	konzentrierte Virus Core-Partikel / concentrated virus core particles
Inaktivierung / Inactivation	Behandlung mit Detergens / treatment with detergent
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:16



Hersteller / Manufacturer:  
Institut Virion\Serion GmbH  
Friedrich-Bergius-Ring 19  
D-97076 Würzburg  
Germany  
Tel. +49 931 3045 0  
Fax +49 931 3045 100  
E-Mail: [dialog@virion-serion.de](mailto:dialog@virion-serion.de)  
Internet: [www.virion-serion.de](http://www.virion-serion.de)