

# Analysenzertifikat

## Certificate of Analysis

### SERION KBR Antigen / CFT Antigen LEGIONELLA PNEUMOPHILA

Best.-Nr. / Order No.	1224
Chargen-Nr. / Lot No.	VDI.AD
Verfallsdatum / Expiry Date	2019-11
Antigen / Antigen	Legionella pneumophila subsp. pneumophila, strain Philadelphia-1
Kultursystem / Culture System	Agar-Plattenkultur / cultivation on agar plates
Präparation / Preparation	beschallter Ganzzellextrakt sonicated cell lysate
Inaktivierung / Inactivation	Zerstörung der Zellen durch Ultraschall / disruption of cells by sonication
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:32



Hersteller / Manufacturer:  
 Institut Virion\Serion GmbH  
 Friedrich-Bergius-Ring 19  
 D-97076 Würzburg  
 Germany  
 Tel. +49 931 3045 0  
 Fax +49 931 3045 100  
 E-Mail: [dialog@virion-serion.de](mailto:dialog@virion-serion.de)  
 Internet: [www.virion-serion.de](http://www.virion-serion.de)