



Analysenzertifikat Certificate of Analysis

virion\sersion

SERION KBR Antigen / CFT Antigen LEPTOSPIRA SEJROE

Best.-Nr. / Order No.	9110
Chargen-Nr. / Lot No.	VHH.AB
Verfallsdatum / Expiry Date	2021-08
Antigen / Antigen	Leptospira interrogans, Serovar sejroe / M84
Kultursystem / Culture System	Flüssigkultur / broth culture
Präparation / Preparation	mit Ultraschall behandeltes Zellmaterial / sonicated cellular material
Inaktivierung / Inactivation	0,1 % Thiomersal / 0.1 % thiomersal
Konservierung / Preservation	0,01 % Thiomersal / 0.01 % thiomersal
Lagerung / Storage	2 °C - 8 °C
Rekonstitution / Reconstitution	1 ml Aqua dest.
Gebrauchsverdünnung / Working Dilution	1:8

serion



multianalyt™

Hersteller / Manufacturer:
Institut Virion\Serion GmbH
Friedrich-Bergius-Ring 19
D-97076 Würzburg
Germany
Tel. +49 931 3045 0
Fax +49 931 3045 100
E-Mail: dialog@virion-serion.de
Internet: www.virion-serion.de